

INDICATION		TAA					
POSITION / LANDMARK		Supine / Chin					
START/END LOCATIONS		Above apex through lungs					
CONTRAST PARAMETERS		50cc, 4cc per sec					
RESPIRATORY PHASE		Inspiration					
SCAN DELAY		SMART PREP tracker @ aortic arch – ascending aorta					
SCAN TYPE		Helical					
KV	mA	Rot Time (sec)	Pitch	Speed (mm/rot)	Noise Index	ASiR	Dose Reduction
120	Smart mA 50-560	0.5	0.984:1	39.37	27		
TECHNIQUE		Using automated exposure control and adjustment of the mA and/or kV according to patient size, radiation dose to be kept as low as reasonably achievable to obtain optimal diagnostic quality images.					
Scans							
Series #	Series	Body Part	DFOV	Thick/Space	Algorithm	Notes	
1	Loc					AP/Lat	
2	Source data	Chest	See note	2x2	STND	Adjust FOV to include bilateral axilla region	
Recons							
Recon source Series #	Recon	Body Part	Thick / Space	Algorithm	W/L	Notes	
2	AX	Chest	5x5	Lung	2000/400	Adjust FOV to include bilateral axilla region	
2	COR	Chest	2x2	STND			
2	SAG	Chest	2x2	STND		Oblique to aorta	

*Please note, recons are displayed as thickness X spacing



CTA CHEST
GE Revolution

Approval: E. Alvarez, MD

rev:1 8/2025

2D / 3D Processing						
Source: Axial with 2x2 standard window						
1. MIP and VR rotation – rotate 360 with 36 images						
Series required in PACS						
Loc, Dose Report, Source data, ALL recons, ALL 3D recons						

ADDITIONAL INSTRUCTIONS:

*Please note, recons are displayed as thickness X spacing