

ADNEXAL OVARIAN WO

Published: 6/2026				Rev: 1				Approval: Dr. Quiroz				Magnet/Vendor: SIEMENS 3.0T		
Seq	Pulse	Dim	Ori	Scan Dir	Phase Enc	FOV	Slc	Gap	TR	TE	Fat Sat	Freq	Phase	Notes
LOCALIZER – FEMALE PELVIS	fl	2D	Sag	R>>L	A>>P	420	6	1.2	7.8	3.69	None	256	205	
AX T2 HASTE MBH	h	2D	Ax	S>>I	A>>P	380	5	1.0	1500	95	None	320	259	
CORONAL T2 HASTE	h	2D	Cor	A>>P	R>>L	360	7	1.4	2000	82	None	320	256	
AXIAL T1 VIBE DIXON IN-OPP	fl	3D	Ax	S>>I	A>>P	360	3	0.6	5.46	1.32/2.55	Dixon	320	224	Send In/Out/Fat/Water
AX T2 FS BH FULL PELVIS	tse	2D	Ax	S>>I	A>>P	380	5	1.5	2500	97	FS	320	224	
SAGITTAL T2	tse	2D	Sag	R>>L	A>>P	230	3	1.0	3420	99	None	320	192	
AXIAL T2	tse	2D	Ax	S>>I	R>>L	230	3	1.0	8250	106	None	320	256	
AXIAL DIFFUSION FB	epse	2D	Ax	S>>I	A>>P	380	5	1.5	5800	56	FS	134	134	
ADC Mapping	derived													Create ADC map and send to PACS