



ANAL FISTULA WWO

Magnet Vendor/Field: SIEMENS 3.0T	Approved By: Dr. Quiroz	Published/Rev: 7/2025 / Rev 2
--	--------------------------------	--------------------------------------

Seq	Pulse	Dim	Ori	Scan Dir	Phase Enc	FOV (mm)	Slc (mm)	Gap (mm)	TR	TE	Fat Sat	Freq	Phase	Notes
localizer	fl	2D	Sag	A>>P	A>>P	420	6.0	1.2	7.8	3.69	None	256	205	
CORONAL T2 HASTE	HASTE	2D	Cor	A>>P	R>>L	360	6.0	1.2	2000	81.0	None	320	256	
SAGITTAL T2 NO FS	tse	2D	Sag	R>>L	A>>P	190	3.0	0.0	4040	94.0	None	256	205	
CORONAL T2 FS	tse	2D	Cor	A>>P	R>>L	190	3.0	0.0	5610	120.0	FS	320	320	
AXIAL T2 NO FS	tseR	2D	Axial	S>>I	R>>L	180	3.0	0.0	7180	98.0	None	320	256	
AXIAL T2 FS	tse	2D	Axial	S>>I	R>>L	180	3.0	0.0	7610	98.0	FS	320	256	
AXIAL T1 VIBE FS PRE	vibe	3D	Axial	S>>I	A>>P	210	3.5	0.7	7.20	3.36	Water Exc.	384	269	
AXIAL T1 VIBE FS POST 1	vibe	3D	Axial	S>>I	A>>P	210	3.5	0.7	7.20	3.36	FS	384	269	create a subtraction (post-pre)
AXIAL T1 VIBE FS POST 2	vibe	3D	Axial	S>>I	A>>P	210	3.5	0.7	7.20	3.36	FS	384	269	create a subtraction (post-pre)
CORONAL T1 VIBE FS POST	vibe	3D	Cor	A>>P	R>>L	210	3.0	0.6	7.20	3.34	FS	384	269	
AXIAL T1 VIBE FS FULL PELVIS STRAIGHT	vibe	3D	Axial	S>>I	A>>P	340	3.5	0.7	6.70	3.05	FS	384	269	

Combine subtracted series into one series and send to PACS