

CHEST WALL WO

Magnet Vendor/Field: SIEMENS MAGNETOM 3.0T XQ					Approved by: Dr Gange					Published/Rev: 6.2026 / 1				
Seq	Pulse	Dim	Ori	Scan Dir	Phase Enc	FOV (mm)	Slc (mm)	Gap (mm)	TR (ms)	TE (ms)	Fat Sat	Freq	Phase	Notes
CORONAL T2 IR	IR	2D	Cor	A>>P	F >> H	280	3.5	0.5	4000.0	69.00	IR	320	288	<i>TI 205 ms</i>
CORONAL T1	TSE	2D	Cor	A>>P	F >> H	280	3.5	0.5	500.0	9.00	None	384	346	
AXIAL T1	TSE	2D	Ax	S>>I	A >> P	280	4.5	0.6	500.0	9.00	None	384	307	
AXIAL T2 IR	IR	2D	Ax	S>>I	A >> P	280	4.5	0.6	4020.0	54.00	IR	384	269	<i>TI 205 ms</i>
SAGITTAL T2 IR	IR	2D	Sag	R>>L	H >> F	280	4.5	0.6	4190.0	61.00	IR	320	240	<i>TI 205 ms</i>