



ENTEROGRAPHY

Magnet Vendor/Field: Siemens MAGNETOM 3.0T XQ					Approved by: Dr Quiroz					Published/Rev: 6.2025 / 1				
Seq	Pulse	Dim	Ori	Scan Dir	Phase Enc	FOV (mm)	Slc (mm)	Gap (mm)	TR (ms)	TE (ms)	Fat Sat	Freq	Phase	Notes
LOCALIZER	VIBE	2D	Sag	R>>L	A>>P	400	6.0	1.2	7.8	3.69	None	256	179	
CORONAL T2 HASTE BH	HASTE	2D	Cor	A>>P	R>>L	380	6.0	1.2	1300	87	None	320	320	
CORONAL T2 TRUFISP	tfi	2D	Cor	A>>P	R>>L	400	7.0	0	3.50	1.54	None	320	256	
CORONAL T1 VIBE PRE	VIBE	3D	Cor	A>>P	F>>H	400	4.0	0.8	2.93	1.16	SPAIR	384	269	
CARE BOLUS SAG FLUORO	VIBE	2D	Sag	R>>L	H>>F	360	34.0	6.8	32.13	1.68	None	320	237	
CORONAL T1 VIBE IMMEDIATE	VIBE	3D	Cor	A>>P	F>>H	400	4.0	0.8	2.93	1.16	SPAIR	384	269	
CORONAL T1 VIBE #2 POST	VIBE	3D	Cor	A>>P	F>>H	400	4.0	0.8	2.93	1.16	SPAIR	384	269	
CORONAL T1 VIBE #3 POST	VIBE	3D	Cor	A>>P	F>>H	400	4.0	0.8	2.93	1.16	SPAIR	384	269	
AXIAL T2 FS BLADE ABD	TSE	2D	Ax	S>>I	A>>P	360	5.0	1.5	2650	95	FS (Strong)	256	256	
AXIAL T2 FS BLADE PEL	TSE	2D	Ax	S>>I	A>>P	360	5.0	1.5	2650	95	FS (Strong)	256	256	
AXIAL T1 VIBE POST ABD	VIBE	3D	Ax	S>>I	A>>P	360	3.0	0.6	3.60	1.40	FS	320	269	
AXIAL T1 VIBE POST PEL	VIBE	3D	Ax	S>>I	A>>P	360	3.0	0.6	3.60	1.40	FS	320	269	