



LIVER WO

Magnet Vendor/Field: Siemens MAGNETOM 3.0T XQ					Approved by: Dr Quiroz					Published/Rev: 4.2026 / 1				
Seq	Pulse	Dim	Ori	Scan Dir	Phase Enc	FOV (mm)	Slc (mm)	Gap (mm)	TR (ms)	TE (ms)	Fat Sat	Freq	Phase	Notes
LOCALIZER	VIBE	2D	Sag	R>>L	A >> P	400	6.0	1.2	9.5	4.36	None	256	179	
CORONAL T2 HASTE MBH	HASTE	2D	Cor	A>>P	R >> L	380	6.0	1.2	1300	87	None	320	320	
AXIAL T2 FS BLADE NAV	TSE	2D	Ax	S>>I	R >> L	380	5.0	1.5	3000	84	FS(Weak)	256	256	
AX T2 HASTE MBH	HASTE	2D	Ax	S>>I	A >> P	380	5.0	1.0	1500	95	None	320	259	
AXIAL T1 VIBE DIXON	VIBE	3D	Ax	S>>I	A >> P	380	3.0	0.6	4.02	1.29 / 2.52	Dixon	320	240	Send In/Out/Fat/Water to PACS
AXIAL T1 VIBE	VIBE	3D	Ax	S>>I	A >> P	380	3.0	0.6	3.60	1.38	FS	320	269	
AXIAL DIFFUSION	Diff	2D	Ax	S>>I	A >> P	380	5.0	1.5	8100	55	FS (Strong)	134	134	b50, b800
ADC MAPPING	ADC Mapping													Create an ADC Map and send to PACS
CORONAL T1 VIBE	VIBE	3D	Cor	A>>P	R >> L	420	3.5	0.7	3.66	1.25	None	320	224	