



LUNG WWO

Magnet Vendor/Field: Siemens MAGNETOM 3.0T XQ					Approved by: Dr Gange					Published/Rev: 6.2025 / 1				
Seq	Pulse	Dim	Ori	Scan Dir	Phase Enc	FOV (mm)	Slc (mm)	Gap (mm)	TR (ms)	TE (ms)	Fat Sat	Freq	Phase	Notes
LOCALIZER	VIBE	2D	Cor	A>>P	R>>L	400	7.0	1.4	7.4	3.01	None	320	240	
LOCALIZER_CHEST	VIBE	2D	Cor	A>>P	R>>L	450	8.0	2.5	7.4	3.01	None	320	240	
CORONAL T2 HASTE BH	HASTE	2D	Cor	A>>P	F>>H	420	8.0	2.4	1800	84	None	256	192	
AXIAL T2 FS BH	TSE	2D	Ax	S>>I	A>>P	400	5.0	1.5	2070	103	FS (Strong)	320	240	
AXIAL T1 IP OP	VIBE	2D	Ax	S>>I	A>>P	400	5.0	1.5	118	2.46 / 6.15	None	320	218	
AXIAL T1 TSE	TSE	2D	Ax	S>>I	A>>P	400	5.0	1.5	500	10	None	256	179	
AXIAL T1 VIBE FS PRE	VIBE	3D	Ax	S>>I	A>>P	400	3.5	0.7	2.91	1.15	FS	384	288	
AXIAL T1 VIBE FS POST	VIBE	3D	Ax	S>>I	A>>P	400	3.5	0.7	2.91	1.15	FS	384	288	Create subtraction (post-pre)
CORONAL T1 VIBE FS DELAYED	VIBE	3D	Cor	A>>P	F>>H	400	4.0	0.8	2.87	1.16	FS	384	269	

Combine subtractions into one series.