

INDICATION		Neck pain						
POSITION / LANDMARK		Supine / SN						
START/END LOCATIONS		EAM through pedicle of T1. Caudal cranial						
CONTRAST PARAMETERS		na						
RESPIRATORY PHASE		na						
SCAN DELAY		na						
SCAN TYPE		Helical						
	KV	mA	Rot Time (sec)	Pitch	Speed (mm/rot)	Noise Index	ASiR	Dose Reduction
	120	Smart mA 50-500	0.8	0.969:1	19.37	16		
TECHNIQUE		Using automated exposure control and adjustment of the mA and/or kV according to patient size, radiation dose to be kept as low as reasonably achievable to obtain optimal diagnostic quality images.						
Scans								
Series #	Series	Body Part	DFOV	Thick/Space	Algorithm	Notes		
1	Loc					AP/Lat		
2	Source data	CSP	25	2 x 2	Bone 3500/800	common angle with spinal cord, bone and STND should be reconstructed together to match		
Recons								
Recon source Series #	Recon	Body Part	Thick / Space	Algorithm	W/L	Notes		
2	AX	CSP	2x2	STND	350/40	common angle with spinal cord, bone and STND should be reconstructed together to match		
2	COR	CSP	2x2	Bone	3500/800			

*Please note, recons are displayed as thickness X spacing

2	COR	CSP	2x2	STND	350/40	
2	SAG	CSP	2x2	Bone	3500/800	
2	SAG	CSP	2x2	STND	350/40	
2D / 3D Processing						
Series required in PACS						
Loc, ALL recons, Dose Report						

ADDITIONAL INSTRUCTIONS: