



Sequence	Pulse Sequence	FOV	Slice Thickness (mm)	Gap (mm)	Matrix		Scan Direction	NOTES
					Phase	Freq		
3 plane loc								
SAG T2	TSE	23	3	0.5	250	384	L-R	No FS
SAG T1	TSE	23	3.0	0.5	268	384	L-R	No FS
SAG STIR	IR	23	3.0	0.5	250	384	L-R	
AX T2 3D	MEDIC	18	3.0	0	256	256	S-I	No FS GE = MERGE PHILLIPS=MFFE HITACHI=ADAGE
CONTRAST								
AX 3D T1	VIBE	18	3.0	0	187	256	S-I	No FS
AX 3D T1 +C	VIBE	18	3.0	0	187	256	S-I	No FS
SAG T1 +C	TSE	23	3.0	0.5	269	384	R-L	FS Strong