

URETHRA DIVERTICULA WWO

Published: 6/8/2026				Rev: 1				Approval: Dr Quiroz			Magnet/Vendor: SIEMENS MAGNETOM 3.0T XQ			
Seq	Pulse	Dim	Ori	Scan Dir	Phase Enc	FOV (mm)	Slc (mm)	Gap (mm)	TR	TE	Fat Sat	Freq	Phase	Notes
Localizer	fl	2D	Sag	R>>L	A>>P	420	6.0	1.2	7.8	3.69	None	256	205	
Coronal T2 HASTE	haste	2D	Cor	A>>P	R>>L	360	6.0	1.2	2000	81	None	320	256	
Axial T2 NO FS	tse	2D	Ax	S>>I	R>>L	190	4.0	0.5	3000	100	None	320	256	
Axial T2 FS HASTE	haste	2D	Ax	S>>I	A>>P	390	5.0	1.5	1600	95	FS	320	256	
Sagittal T2 NO FS	tse	2D	Sag	R>>L	A>>P	190	4.0	0.5	3920	98	None	320	256	
Coronal T2 NO FS	tse	2D	Cor	A>>P	R>>L	190	4.0	0.5	3200	101	None	320	256	
Axial T1 VIBE Dixon	vibe	3D	Ax	S>>I	A>>P	370	3.0	0	5.74	2.46/3.69	Dixon	320	224	Send In/Out/Fat/Water to PACS
Axial T1 FS PRE	vibe	3D	Ax	S>>I	A>>P	370	3.5	0	6.7	2.98	FS	384	268	
Axial T1 FS POST #1	vibe	3D	Ax	S>>I	A>>P	370	3.5	0	6.7	2.98	FS	384	268	create a subtraction (post-pre)
Axial T1 FS POST #2	vibe	3D	Ax	S>>I	A>>P	370	3.5	0	6.7	2.98	FS	384	268	create a subtraction (post-pre)
Axial T1 FS POST #3	vibe	3D	Ax	S>>I	A>>P	370	3.5	0	6.7	2.98	FS	384	268	create a subtraction (post-pre)
Axial T1 FS POST #4	vibe	3D	Ax	S>>I	A>>P	370	3.5	0	6.7	2.98	FS	384	268	create a subtraction (post-pre)

Combine subtractions into one series and sent to PACS